

ELASTY

Cross-linked HA filler

Main Composition: Cross-Linked Hyaluronic Acid, Lidocaine

Country of Origin

REPUBLIC OF KOREA

Product Range



ELASTY D Plus 1cc

1 Syringes of Deep with Lidocaine



ELASTY F Plus 1cc

1 Syringes of Fine with Lidocaine

ELASTY



ELASTY G Plus 1cc

1 Syringes of Grand with Lidocaine

Details about ELASTY

ELASTY: Premium Hyaluronic Acid Dermal Filler for Natural and Long-Lasting Enhancements

ELASTY is a high-quality dermal filler developed by Dongbang Medical Co., Ltd. in South Korea. Utilizing patented Low-Temperature Stabilization Reaction (LTSR) technology, ELASTY fillers are designed to provide exceptional stability, elasticity, and longevity. Widely used for wrinkle correction, volume restoration, and facial contouring, ELASTY ensures natural and precise results.

ELASTY: F PLUS

The softest and most fluid option in the line. It's designed for Superficial or Fine wrinkles and delicate/sensitive areas where you want subtle, natural results without much projection. It integrates well into thinner skin and is less likely to cause lumps in mobile areas.

ELASTY: D PLUS

The middle-ground with medium viscosity. It offers a good balance of elasticity and support, making it suitable for moderate to deep wrinkles and moderate volume restoration.

ELASTY: G PLUS

The thickest and densest variant with the largest particles and highest viscosity. It's built for maximum volume effect and structural support in deeper layers.

ELASTY

Name	Figures	Product specification			Needle specification		
		HA	PHASIC (GEL TYPE)	BDDE	CONTENTS	GAUGE	QUANTITY
ELASTY Grand		24 mg/ml	MONO	<0.2ppm	1ml *1syringe	27G	1EA
ELASTY Deep		24 mg/ml	MONO	<0.2ppm	1ml *1syringe	27G	1EA
ELASTY Fine		24 mg/ml	MONO	<0.2ppm	1ml *1syringe	27G	1EA

Key Features of ELASTY

◦ Primary Ingredients:

Cross-Linked Hyaluronic Acid (HA): Derived from non-animal sources and manufactured under stringent purification processes, HA hydrates, volumizes, and enhances elasticity.

Lidocaine: Included in ELASTY Plus products to minimize discomfort during the procedure.

◦ Effect Duration

Longevity: Results last 6-12 months, depending on the area treated and individual metabolism.

Immediate Results: Enhancements are noticeable right after the procedure.

Safety and Precautions for ELASTY

ELASTY

To ensure safety and optimal results, adhere to these guidelines:

Qualified Practitioner: Ensure the procedure is performed by a licensed and experienced professional.

Medical History: Share any allergies, previous filler treatments, or underlying health conditions with your practitioner.

Contraindications:

- Not recommended during pregnancy or while nursing.
- Avoid treatment if there are active infections, rashes, or inflammation in the target area.

Pre-Treatment Care: Inform your provider about blood-thinning medications or supplements to minimize bruising or bleeding risks.

Post-Treatment Care:

- Avoid touching, rubbing, or applying pressure to the treated area for 24 hours.
- Refrain from strenuous activities, excessive alcohol consumption, and exposure to heat for 48 hours.
- Avoid applying makeup or using harsh skincare products for at least 12-24 hours.

Common Side Effects: Temporary redness, swelling, or bruising at the injection site is normal and typically resolves within a few days.

Maintenance Treatments: Regular follow-up sessions every 6-12 months are recommended to maintain optimal results.

Achieve Radiant Results with ELASTY

ELASTY offers a reliable solution for achieving natural, youthful-looking results with minimal downtime. With its advanced formulation and adherence to safety protocols, ELASTY helps patients enhance their appearance and confidence through precise and long-lasting facial rejuvenation.

For more information visit [ELASTY - PRODUCT](#)