


VOM V is a premium hyaluronic acid-based dermal filler developed by CG Bio Co., Ltd. Designed to address fine lines and initial signs of aging, it utilizes a unique multi-layered phasic structure and R² (Rotation & Revolution) technology. This formulation ensures optimal viscosity and stability, allowing for precise application and long-lasting results. The inclusion of 0.3% lidocaine enhances patient comfort during procedures.

Formulation

- **Hyaluronic Acid (HA):** 20 mg/mL of cross-linked polymeric hyaluronic acid, providing hydration and volume.
- **Lidocaine Hydrochloride:** 0.3% concentration to minimize discomfort during injections.
- **Multi-Layered Phasic Structure:** Combines gel and particle phases for enhanced stability and reduced migration.



- **R² Technology:** Ensures uniform distribution and consistent results through rotational and revolutionary mixing processes.

Benefits

- **Fine Wrinkle Correction:** Effectively smooths superficial facial wrinkles, including crow's feet, perioral lines, and fine lines in the neck area.
- **Skin Hydration:** Enhances moisture retention, leading to improved skin texture and elasticity.
- **Natural-Looking Results:** Provides subtle enhancements that maintain facial harmony.
- **Minimal Downtime:** Quick procedure with immediate visible improvements and minimal recovery time.

Indications

VOM V is indicated for:

- Correction of fine facial wrinkles (e.g., crow's feet, perioral lines).
- Enhancement of lip contour and volume.
- Improvement of skin texture and hydration in the neck area.

Scope of Application

This filler is intended for injection into the intradermal layer of the skin. Recommended treatment areas include:

- Periorbital region (around the eyes)
- Perioral region (around the mouth)
- Neck
- Lips

Administration should be performed by trained medical professionals using appropriate injection techniques.

Duration of Effect

- **Onset:** Immediate visible improvement post-injection.
- **Longevity:** Effects typically last between 6 to 12 months, varying based on individual metabolism, lifestyle, and treatment area.
- **Maintenance:** Follow-up treatments can be scheduled to maintain desired results.

Packaging and Specifications

- **Packaging:** 2 pre-filled syringes containing 1.0 mL each.
- **Needles:** Includes 27G ½" and 30G ½" sterile needles for precise application.
- **Storage:** Store at a temperature between 2°C and 25°C, away from direct sunlight.
- **Shelf Life:** 24 months from the date of manufacture.
- **Manufacturer:** CG Bio Co., Ltd., South Korea.